

# Hydroliq WASSERDESINFEKTION

Disinfectant concentrate for water disinfection



## Water-based disinfectant concentrate

- ~ virucidal, bactericidal, fungicidal
- ~ Dermatest “Excellent”
- ~ Swiss BAG and BAUA-listed
- ~ free of alcohol and synthetic additives
- ~ without caustic or toxic residues
- ~ neutralization of all bacterial odors
- ~ heat and cold resistant
- ~ no requirement of any hazard warnings

## Product Description

Hydroliq WASSERDESINFEKTION is an electrochemically modified water without any synthetic additives. The manufacturing process creates a surface disinfection solution with the ability to eliminate 99.99% of all bacteria, viruses, germs, spores and fungi within seconds.

The chlorine alkali electrolysis process produces a disinfectant solution that is pH-neutral and at the same time very stable (guaranteed 12 months after opening).

## Unique Hygiene Technology

In contrast to conventional room or surface disinfectants, Hydroliq WASSERDESINFEKTION does not leave any residues and thus do not damage any materials during disinfection.

Since no alcohol is used and due to the good skin tolerance the user's skin does not dry out, which prevents germs from penetrating into cracked skin areas.

Hydroliq WASSERDESINFEKTION provides a large spectrum of efficacy and deals with the majority of all germs that are adhering on surfaces:

- ~ **Virucidal**  
HBV, HCV, HIV, influenza, Coronaviridae
- ~ **Bactericidal**  
Staphylococcus aureus, Enterococcus hirae, Pseudomonas aeruginosa, etc.
- ~ **Fungicidal**  
Candida albicans, Aspergillus brasiliensis

## Areas of Application

Hydroliq WASSERDESINFEKTION is used in all hygiene-relevant areas:

- ~ healthcare facilities
- ~ cosmetic- and tattoo studios
- ~ nursing homes
- ~ public buildings
- ~ gastronomy and hotel industry
- ~ industry and commerce
- ~ veterinary sector
- ~ home care

# Hydroliq WASSERDESINFEKTION

Disinfectant concentrate for water disinfection



## Advice of Application

Disinfectant concentrate for disinfecting / treating drinking water.

### Dosage recommendations:

Drinking water e.g. in medical / dental practices: 0.05 – 0.15%

Water in swimming pools and whirlpools: 0.1 – 0.6%

Drinking water in animal husbandry: up to approx. 1.5%

The appropriate dosage is determined taking into account the type and condition of the water installations, as well as the ordinance of the EDI on drinking water and water in public accessible bathrooms and shower facilities. For professional and private use.

Area of application	Dosage	Contact time
Virus inactivation (HBV, HCV, HIV, Influenza, coronaviridae)	undiluted	2 min.
Bacteria inactivation (Staphylococcus aureus, Enterococcus hirae, Pseudomonas aeruginosa, etc.)	undiluted	1 min.
Fungi inactivation (Candida albicans, Aspergillus brasiliensis)	undiluted	1 min.

## Certification

Hydroliq WASSERDESINFEKTION has been dermatologically tested and found to be „EXCELLENT“ by Dermatest.

Hydroliq WASSERDESINFEKTION is Halal certified.

## Packing Unit

Size	Item No.	Content
Hydroliq WASSERDESINFEKTION 5L	HydroW1005	96 pcs/ pallet
Hydroliq WASSERDESINFEKTION 20L	HydroW10020	24 pcs/ pallet
Hydroliq WASSERDESINFEKTION 500L	HydroW100500	
Hydroliq WASSERDESINFEKTION 1000L	HydroW1001000	

## Composition

Active chlorine, released from hypochlorous acid.

## Physical and Chemical Data

The product is not flammable.

## Indicator according to International Directives

The expiration date of this product is printed on the label. Do not use the product after this date!

Use biocides carefully. Always read the label and product information before use. Safety Datasheet will be provided upon request.

BAUA: N-100031

BAG: CHZN6529

## Storage

Keep container tightly closed in a cool, dry, well-ventilated place. Protect from heat. Protect from light.

## Durability

from date of manufacture (expiration date)	24 months
in open condition, but firmly resealed containers or in use with dosing pump	12 months
in a dispenser	12 months