



HYDROLIQ
PURE WATER SOLUTIONS

Hydroliq GREEN

Surface disinfection for agriculture



Water based disinfection

- ~ bactericidal, fungicidal, sporicidal
- ~ Dermatest “Excellent”
- ~ Swiss BAG / BAUA-listed
- ~ free of ethanol and synthetic additives
- ~ without caustic or toxic residues
- ~ heat and cold resistant
- ~ no requirement of any hazard warnings

Product Description

Hydroliq GREEN is an electrochemically modified water without any synthetic additives. The manufacturing process creates a surface disinfection solution with the ability to eliminate 99.99% of all bacteria, germs, spores and fungi within seconds.

The chlorine alkali electrolysis process produces a disinfectant solution that is pH-neutral and at the same time very stable (guaranteed 12 months after opening).

Unique Hygiene Technology

In contrast to conventional disinfectants, Hydroliq GREEN does not leave any caustic residues and therefore does not damage any materials during disinfection.

Since no alcohol is used and due to the good skin tolerance the user's skin does not dry out, which prevents germs from penetrating into cracked skin areas.

Hydroliq GREEN provides a large spectrum of efficacy and deals with the majority of all germs that are adhering on surfaces:

Areas of Application

Hydroliq GREEN is used in all agricultural areas:

- ~ Feed storage
- ~ Feed sector
- ~ greenhouses
- ~ foil tunnel
- ~ farms
- ~ horticulture
- ~ sports venues
- ~ cutting devices



HYDROLIQ
PURE WATER SOLUTIONS

Hydroliq GREEN

Surface disinfection for agriculture



Advice of Application

Disinfection of non-porous surfaces indoors and outdoors for agriculture, Greenhouses, feed storage, the food sector and for the private sector. Application by wiping or spraying. Before disinfection is a cleaning of the surfaces required. Use a sufficient amount so that the treated areas remain moist during the entire duration of exposure. For professional and private use.

Area of application	Dosage	Contact time
Bacteria inactivation	undiluted	1 min.
Fungi / Spores inactivation	undiluted	1 min.

Certification

Hydroliq GREEN *Concentrate* has been dermatologically tested and found to be „EXCELLENT“ by dermatest.

Packing Unit

Size	Item No.	Content
Hydroliq GREEN 5L	HydroGr1005	96 pcs./ pallet
Hydroliq GREEN 20L	HydroGr10020	28pcs./ pallet
Hydroliq GREEN 500L	HydroGr100500	
Hydroliq GREEN 1000L	HydroGr1001000C	

Composition

Active chlorine, released from hypochlorous acid

Physical and Chemical Data

The product is not flammable.

Indicator according to International Directives

The expiration date of this product is printed on the label. Do not use the product after this date!

Use biocides carefully. Always read the label and product information before use. Safety Datasheet will be provided upon request.

BAUA: N-96237

BAG: CHZN6380

Storage

Keep container tightly closed in a cool, dry, well-ventilated place. Protect from heat. Protect from light.

Durability

The following durability data should be observed when using Hydroliq GREEN:

from date of manufacture (expiration date)	24 months
in open condition, but firmly resealed containers or in use with dosing pump	12 months