

# Hydroliq ANIMAL Concentrate

Disinfection for veterinary section



**HYDROLIQ**  
PURE WATER SOLUTIONS



## Water based disinfection

- ~ Virucidal, bactericidal, fungicidal
- ~ Dermatest “Excellent”
- ~ Swiss BAG / BAUA-listed
- ~ free of ethanol
- ~ without caustic residues
- ~ no requirement of any hazard warnings
- ~ hardly forms any resistance

## Product Description

Hydroliq ANIMAL Concentrate is an electrochemically modified water. Our manufacturing process creates a surface disinfection solution with the ability to eliminate 99.99% of all virus, bacteria, germs, spores and fungi.

The chlorine alkali electrolysis process produces a disinfectant solution that is pH-neutral and at the same time very stable (guaranteed 12 months after opening).

## Unique Hygiene Technology

In contrast to conventional disinfectants in the veterinary sector, the solutions from Hydroliq Pure Water Solutions are non-corrosive and do not leave residues.

Since no alcohol is used and due to the good skin tolerance the user's skin does not dry out, which prevents germs from penetrating into cracked skin areas.

Hydroliq ANIMAL Concentrate provides a large spectrum of efficacy and deals with the majority of all germs that are adhering on surfaces:

- ~ **Virucidal**  
HBV, HCV, HIV, Influenza, Coronaviridae, Influenza-A viruses (Bird flu), Asfarviridae, Genus Asfivirus ASPV-Virus (Swine fever), etc.
- ~ **Bactericidal**  
Staphylococcus aureus, Enterococcus hirae, Pseudomonas aeruginosa, Dichelobacter nodosus, Fusobacterium necrophorum, Streptococcus agalactiae, Streptococcus dysgalactiae, Streptococcus uberis, Mycoplasmosis, E. coli
- ~ **Fungicidal**  
Candida albicans, Aspergillus brasiliensis

## Areas of Application

Hydroliq ANIMAL Concentrate is used in all veterinary areas:

- ~ chicken farming
- ~ animal care
- ~ mast operation
- ~ rearing

# Hydroliq ANIMAL Concentrate

Disinfection for veterinary section



## Advice of Application

For the preventive disinfection of empty cages, stalls and terrariums. We recommend a diluted application (1:50 to 1:15) for Spray, wipe or dip processes. For commercial and private use. Eliminates fungi, spores and all kinds of germs.

## Certification

Hydroliq ANIMAL *Concentrate* has been dermatologically tested and found to be „EXCELLENT“ by dermatest.

## Packing Unit

Size	Item No.	Content
Hydroliq ANIMAL <i>Concentrate</i> 5L	HydroAn1005	96 pcs./ pallet
Hydroliq ANIMAL <i>Concentrate</i> 20L	HydroAn10020	28 pcs./ pallet
Hydroliq ANIMAL <i>Concentrate</i> 500L	HydroAn100500C	

## Composition

Active chlorine, released from hypochlorous acid

## Physical and Chemical Data

The product is not flammable.

## Indicator according to International Directives

The expiration date of this product is printed on the label. Do not use the product after this date!

Use biocides carefully. Always read the label and product information before use. Safety Datasheet will be provided upon request.

BAUA: N-96234

BAG: CHZN5775.03.001

## Storage

Keep container tightly closed in a cool, dry, well-ventilated place. Protect from heat. Protect from light.

## Durability

The following durability data should be observed when using Hydroliq ANIMAL *Concentrate*:

from date of manufacture  
(expiration date) 24 months

in open condition, but firmly  
resealed containers or in use with  
dosing pump 12 months