



HYDROLIQ
PURE WATER SOLUTIONS

Hydroliq GREEN *Ready to use*

Surface disinfection for agriculture



Water based disinfection

- ~ bactericidal, fungicidal, sporicidal
- ~ Dermatest "Excellent"
- ~ Swiss BAG / BAUA-listed
- ~ free of ethanol
- ~ without caustic residues
- ~ no requirement of any hazard warnings

Product Description

Hydroliq GREEN *Ready to use* is an electrochemically modified water. Our manufacturing process creates a surface disinfection solution with the ability to eliminate 99.99% of all bacteria, germs, spores and fungi.

The chlorine alkali electrolysis process produces a disinfectant solution that is pH-neutral and at the same time very stable (guaranteed 12 months after opening).

Unique Hygiene Technology

In contrast to conventional disinfectants, Hydroliq GREEN *Ready to use* does not leave any caustic residues and therefore does not damage any materials during disinfection.

Since no alcohol is used and due to the good skin tolerance the user's skin does not dry out, which prevents germs from penetrating into cracked skin areas.

Hydroliq GREEN *Ready to use* provides a large spectrum of efficacy and deals with the majority of all germs that are adhering on surfaces:

- ~ **Bactericidal**
- ~ **Viruzidal**
- ~ **Fungicidal**
- ~ **Sporicidal**

Areas of Application

Hydroliq GREEN *Ready to use* is used in all agricultural areas for non-planted surfaces:

- ~ greenhouses
- ~ foil tunnel
- ~ farms
- ~ horticulture
- ~ sports venues
- ~ cutting devices



HYDROLIQ
PURE WATER SOLUTIONS

Hydroliq GREEN *Ready to use*

Surface disinfection for agriculture



Advice of Application

Use for the surface disinfection in sports facilities with floors of sand, gravel, bitumen, tar and synthetic materials, disinfection of synthetic turf, footways as well as parts of buildings and cutting devices. The undiluted product can be applied by spraying, wiping or immersion (up to three times per week recommended). For professional and private use. Eliminates fungi, spores and all types of germs.

Area of application	Dosage	Contact time
Bacteria inactivation	undiluted	1 min.
Fungi and Spores inactivation	undiluted	1 min.

Certification

Hydroliq GREEN *Ready to use* has been dermatologically tested and found to be „EXCELLENT“ by dermatest.

Packing Unit

Size	Item No.	Content
Hydroliq GREEN <i>Ready to use</i> 1L	HydroGr0101	12 pcs./ box
Hydroliq GREEN <i>Ready to use</i> 5L	HydroGr0105	96 pcs./ pallet

Composition

Active chlorine, released from hypochlorous acid

Physical and Chemical Data

The product is not flammable.

Indicator according to International Directives

The expiration date of this product is printed on the label. Do not use the product after this date!

Use biocides carefully. Always read the label and product information before use. Safety Datasheet will be provided upon request.

BAUA: N-96238

BAG: CHZN5775.02.002

Storage

Keep container tightly closed in a cool, dry, well-ventilated place. Protect from heat. Protect from light.

Durability

The following durability data should be observed when using Hydroliq GREEN *Ready to use*:

from date of manufacture (expiration date)	24 months
in open condition, but firmly resealed containers or in use with dosing pump	12 months