

Hydroliq AIR Pro

Disinfection by fogging



HYDROLIQ
PURE WATER SOLUTIONS

Disinfection by fogging

- › Spectrum of activity: virucidal, bactericidal, fungicidal, sporicidal
- › Listed by BAG / ECHA / VAH
- › Free from harmful substances and alcohol
- › Gentle on materials
- › Odour-neutralizing
- › No labeling required according to GHS



Product Description

Discover the uncomplicated method of disinfection through nebulization! Fogging offers a simple and effective solution for disinfecting indoor spaces. From operating theatres, patient rooms and surgeries to busy classrooms and concert halls – with our nebulization technology you can create a controlled and safe indoor climate everywhere.

Experience the simplicity of fogging disinfection without compromising on effectiveness.

Application

Airtransmitted room disinfection by cold fogging. Application is carried out using a device without the presence of people in the room. After the exposure time: Wait 1 minute and ventilate well before entering the treated room.

Area of application

Operating theater
Surgery rooms
Patient rooms
Classrooms

Storage

Store the product tightly closed in a cool, dry place.
Protect from light and heat.

*Handle biocides with care. Always read the label and product information before use.
Safety data sheet available on request.*

Certification

Approval number: CHZN5775.04.002



295687 - CH - 0202

Composition

Active chlorine released from hypochlorous acid.
1000ppm

Effect

Bactericidal (EN13727)
Myko Bactericidal (EN14348)
Virucidal (EN14476)
Fungicidal (EN13624)
Sporicidal (EN17126)

Hazard warnings

No instructions necessary

Packaging Units

Size	Item number	Content
Hydroliq AIR Pro 5L	HydroAi5005K	96 pieces/pallet
Hydroliq AIR Pro 20L	HydroAi100020K	24 pieces/pallet
Hydroliq AIR Pro 500L	HydroAi1000500C	
Hydroliq AIR Pro 1000L	HydroAi10001000C	